

BIOSINKER II

By HPNE

Gamma-stable PVDF Tube Weight



For the applications requiring the highest purity and strength, our Biosinker tube weights the ideal solution. Constructed of PVDF, our Biosinker II will maintain optimal depth and placement of your tube to obtain the most efficient flow possible.

Biosinker II represents a slight design change from Biosinker I, featuring a new tapered bottom allowing it to reach every space in a bottle while maintaining optimal flow. The Biosinker II is autoclave and gamma-stable, but can also be provided non-irradiated.

Uses include:

- A direct replacement for a rigid dip tube
- Keeps tube weighted during the dispense cycle
- Prevents curling, sticking or floating

PART #	HOS BARB	DIAMETER
HPN01-BSKN-125	1/8"	0.063"
HPN01-BSKN-250	1/4"	0.188"
HPN01-BSKN-375	3/8"	0.313"
HPN01-BSKN-500	1/2"	0.438"



hp-ne.com | info@hp-ne.com | 401.349.4477

High Purity New England Inc., 2 Thurber Boulevard, Smithfield, RI 02917